

## **Laboratory Medicine - Test Information Guide Chlamydia trachomatis Culture**

Alternate Name:

Domontosonti	Missahialass
Department:	Microbiology
STAT Test Available:	N
TAT (Days):	
Specimen Type:	MultiTrans Media
Collection Information & Patient Prep:	
<ul> <li>Cultures performed for test-of- cure or alleged/suspected sexual assault. Do not use expired medium or medium which has turned yellow or purple or which appears to contain particulate matter.</li> <li>Procedure:         <ol> <li>Reasonable vigorous swabbing is necessary to collect cellular material, since viable organisms are usually intracellular.</li> <li>For cervical specimens, discharge should be wiped away before collecting the sample. The discharge is unlikely to contain viable organisms and may be toxic to the culture. The squamocolumnar epithelial junction is the area of the female genital tract most likely to yield the agent.</li> </ol> </li> </ul>	
Methodology:	
Requisition:	PHL
Critical Result Ranges:	
Testing Referred To:	PHL
Sample Volume:	
Unit of Measure:	
Image:	
Reference Ranges:	