

## Laboratory Medicine - Test Information Guide

### Tetanus Serology

Alternate Name:

Department: Microbiology

STAT Test Available: N

TAT (Days):

Specimen Type: Serum separator tube

Collection Information & Patient Prep:

- Testing is available only for the rare event of an adverse reaction to the Tetanus vaccine or the possibility of humoral immunodeficiency in the patient. This must be indicated on the test requisition in order for testing to be performed.
- Store in 4oC refrigerator for up to 48 hours.

Methodology:

Requisition: PHL General Test Requisition Form

Critical Result Ranges:

Testing Referred To: PHL

Sample Volume: 1.0 mL serum

Unit of Measure:

Image:

Reference Ranges: