

Laboratory Medicine - Test Information Guide

Toxoplasmosis serology

Alternate Name: Neonatal and Adult testing

Department: Microbiology

STAT Test Available: N

TAT (Days):

Specimen Type: Serum separator tube

Collection Information & Patient Prep:

- If congenital infection is suspected, submit maternal serum as well as one from the infant.
- Indicate on the requisition if the test is for diagnosis or immunity.
- Store in 4oC refrigerator.

Methodology:

Requisition: PHL

Critical Result Ranges:

Testing Referred To: PHL

Sample Volume: 1.0 mL

Unit of Measure:

Image:

Reference Ranges: