

Laboratory Medicine - Test Information Guide

Urine or Urine Cystoscopy for AFB, TB, MAI

Alternate Name: Urine or Urine Cystoscopy for Mycobacteria

Department: Microbiology

STAT Test Available: N

TAT (Days):

Specimen Type: Sterile screw cap container

Collection Information & Patient Prep:

- Refer to [SCI-5.7.18.1.0](#)
- Catheter or midstream urine: at least 40.0 ml of the first morning specimen.
- Suprapubic tap: at least 40 ml. Submit in a sterile container
- Store in 4oC refrigerator

Methodology:

Requisition: PHL

Critical Result Ranges:

Testing Referred To: PHL

Sample Volume:

Unit of Measure:

Image:

Reference Ranges: