

Laboratory Medicine - Test Information Guide

1p19q CODEL

Alternate Name: N/A

Department: Molecular Diagnostics

STAT Test Available: N/A

TAT (Days): Refer to Reference Laboratory

Specimen Type: 7mL in Sodium Heparin* (GREEN)

Collection Information & Patient Prep:

- TESTING PERFORMED ON SUNNYBROOK PATIENTS ONLY
- Sample must be received Monday - Friday by 10am (excluding STAT holidays) for same day delivery within the GTA, to the Molecular Services Laboratory (E410) for referral to testing facility.
- Transport all specimens in a sealed biohazard bag maintaining a barrier between the specimen and documentation.
- Transport all specimens at ambient temperature.

Methodology: NA

Requisition: [UHN Requisition](#)

Critical Result Ranges:

Testing Referred To:

Sample Volume:

Unit of Measure:

Image:

Reference Ranges: