

Laboratory Medicine - Test Information Guide

BCR-ABL1 Diagnostic

Alternate Name:

Department: Molecular Diagnostics

STAT Test Available: N/A

TAT (Days): Refer to Reference Laboratory

Specimen Type: Peripheral blood - 4 EDTA (Lavender)

Collection Information & Patient Prep:

- TESTING PERFORMED ON SUNNYBROOK PATIENTS ONLY
- Sample must be received Monday - Friday by 10am (excluding STAT Holidays) for same day delivery within the GTA, to the Molecular Services Laboratory (E410) for referral to testing facility.
- Transport all specimens in a sealed biohazard bag maintaining a barrier between the specimen and documentation.
- Transport all specimens at ambient temperature

Methodology:

Requisition: [UHN Requisition](#)

Critical Result Ranges:

Testing Referred To:

Sample Volume: 7ml

Unit of Measure:

Image:

Reference Ranges: