

Laboratory Medicine - Test Information Guide

Biopsies-Tissues for Virus Isolation

Alternate Name: NA

Department: Microbiology

STAT Test Available: N

TAT (Days):

Specimen Type:

Collection Information & Patient Prep:

- Sterile screw cap container or MultiTrans Media
- MultiTrans Media obtained from Microbiology lab.
- Place in virus transport to prevent drying of tissue.
- Indicate the virus queried and diagnosis.
- Refer to [SCI-5.7.15.2.10](#)

Methodology:

Requisition: PHL

Critical Result Ranges:

Testing Referred To: PHL

Sample Volume:

Unit of Measure:

Image:

Reference Ranges: