

Laboratory Medicine - Test Information Guide

Chlamydomphila pneumoniae-Mycoplasma pneumoniae- Respiratory PCR

Alternate Name:

Department: Microbiology

STAT Test Available: N

TAT (Days):

Specimen Type:

Collection Information & Patient Prep:

- 1 - 2 mL fluid from BAL in a sterile screw cap container or nasopharyngeal swab in Multitrans Media
- Specimens should be collected during the acute phase of illness.
- MultiTrans Media obtained from Specimen Management Centre.
- Store in 4oC refrigerator.

Methodology:

Requisition: PHL General Test Requisition Form

Critical Result Ranges:

Testing Referred To: PHL

Sample Volume:

Unit of Measure:

Image:

Reference Ranges: