

Laboratory Medicine - Test Information Guide CSF for Syphilis

Alternate Name:	
Department:	Microbiology
STAT Test Available:	Y Refer to <u>SCI-5.7.15.1.0</u>
TAT (Days):	
Specimen Type:	CSF collected by LP, ventriculostomy, or ommaya fluid in a sterile screw top container
Collection Information & Patient Prep:	Store in 4C refrigerator
Methodology:	
Requisition:	PHL
Critical Result Ranges:	
Testing Referred To:	PHL
Sample Volume:	1 mL
Unit of Measure:	
Image:	
Reference Ranges:	