



Laboratory Medicine - Test Information Guide

Donath-Landsteiner Test

Alternate Name:

Department: Blood and Tissue Bank

STAT Test Available: N

TAT (Days): Available Monday to Friday 8-4 only

Specimen Type: Plain red top plastic vacutainer

Collection Information & Patient Prep: Sample must be kept at 35-37 C immediately after collection

Methodology:

Requisition: Blood and Tissue Bank Order Form

Critical Result Ranges:

Testing Referred To:

Sample Volume: 7 - 10 mL

Unit of Measure:

Image:

Reference Ranges: