

Laboratory Medicine - Test Information Guide

FISH - CLL Panel

Alternate Name:

Department: Molecular Diagnostics

STAT Test Available: N/A

TAT (Days): Refer to Reference Laboratory

Specimen Type: 7mL in Sodium Heparin* (GREEN)

Collection Information & Patient Prep:

- TESTING PERFORMED ON SUNNYBROOK PATIENTS ONLY
- Sample must be received Monday - Friday by 10am (excluding STAT Holidays) for same day delivery within the GTA, to the Molecular Services Laboratory (E410) for referral to testing facility.
- Transport all specimens in a sealed biohazard bag maintaining a barrier between the specimen and documentation.
- Transport all specimens at ambient temperature.

Methodology: Contact Reference Laboratory for Methodology 416-340-4800 X 5725

Requisition: [UHN Requisition](#)

Critical Result Ranges: N/A

Testing Referred To: University Health Network - Cancer Cytogenetics Laboratory

Sample Volume: 7mL

Unit of Measure: Interpretative Comment

Image:

Reference Ranges: