

## Laboratory Medicine - Test Information Guide Hepatitis B HBc IGM

Alternate Name:	
Department:	Microbiology
STAT Test Available:	N
TAT (Days):	
Specimen Type:	Serum separator tube
Collection Information & Patient Prep:	
•	This test is processed only when HBsAg is positive Refer to $\frac{\text{SCI-}5.7.15.2.9}{\text{Store in 4C refrigerator}}$
Methodology:	
Requisition:	PHL
Critical Result Ranges:	
Testing Referred To:	PHL
Sample Volume:	2 mL
Unit of Measure:	
Image:	

Reference Ranges: