

Laboratory Medicine - Test Information Guide

Kidney washing (renal pelvis, pelvic)

Alternate Name: Renal Pelvis washings, Pelvic washing

Department: Cytology

STAT Test Available: Prior Consultation required [SCI-5.7.20.5](#) Anatomic Pathology -STAT Specimen Procedure

TAT (Days): Refer to [SCI-5.7.20.7](#) Anatomic Pathology - Turn Around Times

Specimen Type:

Collection Information & Patient Prep:

- Clinical information is crucial to an accurate diagnosis.
- Use a 90ml sterile collection container.
- Submit fresh (no fixative)
- Specimen must be placed in 4C refrigerator after hours.
- [SCI-5.7.20.8](#) Anatomic Pathology - Specimen identification and requisition
- [SCI-5.7.20.6](#) Anatomic Pathology - Where to deliver the specimen

Methodology:

Requisition:

Critical Result Ranges:

Testing Referred To:

Sample Volume:

Unit of Measure:

Image:

Reference Ranges: