

Laboratory Medicine - Test Information Guide

Synovial Joint Fluid for Virus Isolation

Alternate Name:

Department: Microbiology

STAT Test Available: N

TAT (Days):

Specimen Type: Sterile screw cap container

Collection Information & Patient Prep:

- Indicate the virus queried and diagnosis.
- Store in 4oC refrigerator

Methodology:

Requisition: PHL

Critical Result Ranges:

Testing Referred To: PHL

Sample Volume:

Unit of Measure:

Image:

Reference Ranges: