

Laboratory Medicine - Test Information Guide

Throat swab for Virus Isolation

Alternate Name:

Department: Microbiology

STAT Test Available: N

TAT (Days):

Specimen Type: MultiTrans Media

Collection Information & Patient Prep:

- MultiTrans Media obtained from Specimen Management Centre.
- Indicate the virus queried and diagnosis
- Store in 4°C fridge
- For non-respiratory viruses

Methodology:

Requisition: PHL

Critical Result Ranges:

Testing Referred To: PHL

Sample Volume:

Unit of Measure:

Image:

Reference Ranges: