

Laboratory Medicine - Test Information Guide **BCR-ABL1 MRD**

| Department: | Molecular Diagnostics |
|--|---|
| STAT Test Available: | N/A |
| TAT (Days): | Refer to Reference Laboratory |
| Specimen Type: | EDTA whole blood (Lavender), DO NOT freeze or centrifuge |
| Collection Information & Patient Prep: | |
| • • • | TESTING PERFORMED ON SUNNYBROOK PATIENTS ONLY Sample must be sent to Molecular Diagnostics Laboratory as soon as possible. Sample collected on Friday or before Statutory holidays must be delivered to E410 before noon. Transport all specimens in a sealed biohazard bag maintaining a barrier betwee the specimen and documentation. Transport all specimens at ambient temperature |
| Methodology: | |
| Requisition: | SB Molecular/Biomarkers Requisition |
| Critical Result Ranges: | |
| Testing Referred To: | |
| Sample Volume: | 2 tubes of 4mL |
| Unit of Measure: | |
| Image: | |
| Reference Ranges: | |

barrier between

Alternate Name: